



# CALCIUM

4 x 25 mL, 4 x 50 mL  
51005001, 51005002

## INTENDED USE

This reagent is intended for *in vitro* quantitative determination of calcium in serum, plasma & urine.

- Modified OCPC methodology
- Linear up to 15 mg/dL

## CLINICAL SIGNIFICANCE

Calcium is an important ion present in the body. Mainly it is found in bones. In serum calcium exists equally in a free ionized form & also in a bound form with albumin. Calcium helps in enzyme activation, muscle contraction, coagulation of blood, regulation of some hormonal secretions & cell membrane permeability.

Increased levels are found in hyperthyroidism, malignant tumors, acute osteoporosis & adrenal insufficiency.

Decreased levels found in hypoparathyroidism, osteomalacia, rickets, renal failure & tetanus.

## PRINCIPLE

Calcium OCPC procedure is based on the reaction of calcium ions ( $\text{Ca}^{2+}$ ) with O-cresolphthalein complex in an alkaline solution to form an intense violet coloured complex which shows maximum absorbance at 578nm. The 8-hydroxyquinoline prevents  $\text{Mg}^{2+}$  interference upto 4 mmol/L.

## REAGENT COMPOSITION

<b>CALCIUM DYE REAGENT (R2)</b>	<b>2 x 25 mL / 2 x 50 mL</b>
Diethylamine	360 mmol/L
<b>CALCIUM BASE REAGENT (R1)</b>	<b>2 x 25 mL / 2 x 50 mL</b>
O-Cresolphthalein complex	0.15 mmol/L
8-Hydroxyquinoline	17.2 mmol/L
<b>CALCIUM STANDARD</b>	<b>1 x 4 mL</b>
Calcium standard concentration	10 mg/dL

## STORAGE & STABILITY

The sealed reagents are stable up to the expiry date stated on the label, when stored at 2 - 8°C.

## LINEARITY

This reagent is linear up to 15 mg/dL.

If the concentration is greater than linearity (15 mg/dL), dilute the sample with normal saline and repeat the assay. Multiply the result with dilution factor.

## NORMAL RANGE

It is recommended that each laboratory establish its own reference values.

The following value may be used as guide line.

Serum	: 8.8 - 10.2 mg/dL
Urine	: 100 - 400 mg/24 hrs

## PREPARATION AND STABILITY OF REAGENT

Mix reagent 1 (R1) and Reagent 2 (R2) in the ratio 1:1.

## PRECAUTION

To avoid contamination, use clean laboratory wares.

Avoid direct exposure of reagent to light.

(Use acid washed (50 %  $\text{HNO}_3$ ) glass wares & tips)

## SAMPLE

Serum / plasma (free of haemolysis) / Urine (1/3 diluted)

## GENERAL SYSTEM PARAMETER

Mode of Reaction	End point
Slope of reaction	Increasing
Wavelength	578 nm (565-580 nm)
Temperature	30°C
Standard Concentration	10 mg/dL
Linearity	15 mg/dL
Blank	Reagent
Incubation time	5 min
Sample volume	10 $\mu\text{L}$
Reagent volume	1000 $\mu\text{L}$
Cuvette	1 cm light path

## LABORATORY PROCEDURE

	Blank	Standard	Sample
Working Reagent	1000 $\mu\text{L}$	1000 $\mu\text{L}$	1000 $\mu\text{L}$
Standard	-	10 $\mu\text{L}$	-
Sample	-	-	10 $\mu\text{L}$

Mix and incubate for 5 min. at room temperature. Read the absorbance of standard and sample against reagent blank.

## CALCULATION

$$\text{Calcium Conc. (mg/dL)} = \frac{\text{Absorbance of sample}}{\text{Absorbance of Standard}} \times 10$$

## INTERFERENCE

Bilirubin concentrations higher than 20 mg/dL and phosphate higher than 40 mg/dL, will interfere with the assay.

## BIBLIOGRAPHY

1. Schwarzenbach, G.; *Analyst.*, 80, (1955) 713-729
2. Kessler, G., Wolfman, M., *Clin.Chem.*, 10, (1964) 686 – 703
3. Connerty, H. V., Briggs, A.R., *Am. J. Clin.Pathol.*, 45, (1965) 290-296
4. Gitelmann, H. J. *Anal Biochem* 18, (1967) 521-531
5. Biggs, H.G., Moorehead, W. R.; *Clin.Chem.*, 20, (1974) 1458-1460

SYMBOLS USED ON THE LABELS: IN VITRO DIAGNOSTIC USE SEE PACKAGE INSERT FOR PROCEDURE LOT NUMBER MANUFACTURER'S ADDRESS MANUFACTURING DATE EXPIRY DATE TEMPERATURE LIMIT

**BIOSS BioTechnology, GmbH**  
Boekhulter Weg 1a, 47638 Straelen, Germany  
E-mail: sales@biossbiotech.de, support@biossbiotech.de  
 Arthrex GmbH  
Erwin-Hiescher-Strasse 9, 81249 Munchen, Germany

ADS/V.02/110314

ISO 9001 : 2008  
ISO 13485 : 2003